

# MYOD1 antibody



Catalog Number: 107316

---

## Product name

MYOD1 antibody

## Specificity

Human; other species not tested.

## Antibody description

MYOD1 Mouse Monoclonal antibody. Positive WB detected in human skeletal muscle tissue. Observed molecular weight by Western-blot: 47 kDa

## Preparation

This antibody was obtained by immunization of MYOD1 recombinant protein (Accession Number: NM\_002478). Purification method: Protein G purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



human skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:107316(MYOD1 Antibody) at dilution of 1:1000