

NENF antibody



Catalog Number: 107327

Product name

NENF antibody

Immunogen

[Human NENF Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse; other species not tested.

Antibody description

NENF Mouse Monoclonal antibody. Positive IHC detected in human cerebellum tissue, human colon tissue, human heart tissue, human kidney tissue. Positive WB detected in human brain tissue, human heart tissue. Observed molecular weight by Western-blot: 18 kDa

Preparation

This antibody was obtained by immunization of NENF recombinant protein (Accession Number: NM_013349). Purification method: Protein A purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2a

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:107327(NENF antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human cerebellum using Catalog No:107327(NENF antibody) at dilution of 1:50 (under 10x lens)



NENF antibody



Catalog Number: 107327

Immunohistochemical of paraffin-embedded human cerebellum using Catalog

No:107327(NENF antibody) at dilution of 1:50 (under 40x lens)