

NESP55,GNAS antibody



Catalog Number: 107328

Product name

NESP55,GNAS antibody

Specificity

Human, Pig; other species not tested.

Antibody description

NESP55,GNAS Mouse Monoclonal antibody. Positive WB detected in fetal human brain tissue, pig brain tissue. Observed molecular weight by Western-blot: 46 kDa

Preparation

This antibody was obtained by immunization of NESP55,GNAS recombinant protein (Accession Number: NM_000516). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2b

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



fetal human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:107328(NESP55,GNAS Antibody) at dilution of 1:2000