

NUDT21 antibody



Catalog Number: 107341

Product name

NUDT21 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

NUDT21 Mouse Monoclonal antibody. Positive WB detected in NIH/3T3 cells, HEK-293 cells, HeLa cells, MCF-7 cells. Observed molecular weight by Western-blot: 27 kDa

Preparation

This antibody was obtained by immunization of NUDT21 recombinant protein (Accession Number: NM_007006). Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with (NUDT21 Antibody) at dilution of 1:3000