

# NUMB antibody

Catalog Number: 107342

## Product name

NUMB antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

NUMB Mouse Monoclonal antibody. Positive IF detected in A431 cells. Positive WB detected in Human brain. Observed molecular weight by Western-blot: 66-72 kDa

## Preparation

This antibody was obtained by immunization of peptide (Accession Number: NM\_001005745). Purification method: Caprylic acid/ammonium sulfate precipitation.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2a

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IF: 1:20-1:200

## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:107342(NUMB antibody) at dilution of 1:300



Immunofluorescent analysis of A431 cells using Catalog No:107342(NUMB Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Mouse IgG(H+L)