

# Insulin antibody



Catalog Number: 107350

## Product name

Insulin antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

Insulin Mouse Monoclonal antibody. Positive IHC detected in human pancreas tissue, mouse pancreas tissue.

## Preparation

This antibody was obtained by immunization of Insulin recombinant protein (Accession Number: NM\_001185097). Purification method: Protein A purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2a

## Applications


ELISA, IHC

## Dilutions


Recommended Dilution:

IHC: 1:500-1:5000

## Validations



Immunohistochemistry of paraffin-embedded human pancreas tissue slide using Catalog No:107350(Insulin Antibody) at dilution of 1:5000 (under 10x lens).



Immunohistochemistry of paraffin-embedded human pancreas tissue slide using Catalog No:107350(Insulin Antibody) at dilution of 1:5000 (under 40x lens).