

Integrin beta-1 antibody



Catalog Number: 107352

Product name

Integrin beta-1 antibody

Specificity

Human; other species not tested.

Antibody description

Integrin beta-1 Mouse Monoclonal antibody.
Positive WB detected in HT-1080 cells, HUVEC cells. Observed molecular weight by Western-blot: 130 kDa

Preparation

This antibody was obtained by immunization of Integrin beta-1 recombinant protein (Accession Number: NM_002211). Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:2000-1:20000

Validations



HT-1080 cells were subjected to SDS PAGE followed by western blot with Catalog No:107352(Integrin beta-1 Antibody) at dilution of 1:8000