

# JUN antibody

Catalog Number: 107358

---

## Product name

Monoclonal

JUN antibody

## Ig Type

## Specificity

Mouse IgG1

Human, Mouse; other species not tested.

## Applications

## Antibody description

ELISA, WB

JUN Mouse Monoclonal antibody. Positive WB detected in RAW 264.7 cells, NIH/3T3 cells. Observed molecular weight by Western-blot: 39 kDa

## Dilutions


Recommended Dilution:

WB: 1:1000-1:10000

## Preparation

This antibody was obtained by immunization of JUN recombinant protein (Accession Number: NM\_002228). Purification method: Protein G purified.

## Validations



## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

RAW 264.7 cells were subjected to SDS PAGE followed by western blot with Catalog No:107358(JUN Antibody) at dilution of 1:4000

## Clonality