

# LDHA antibody



Catalog Number: 107375

## Product name

LDHA antibody

## Specificity

Human; other species not tested.

## Antibody description

LDHA Mouse Monoclonal antibody. Positive WB detected in pig liver tissue. Observed molecular weight by Western-blot: 28 kDa, 40 kDa

## Preparation

This antibody was obtained by immunization of LDHA recombinant protein (Accession Number: NM\_005566). Purification method: Protein A purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2b

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



pig liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:107375(LDHA Antibody) at dilution of 1:1500