

# IL-17F antibody



Catalog Number: 107385

## Product name

IL-17F antibody

## Immunogen

## Specificity

Human, Mouse; other species not tested.

## Antibody description

IL-17F Mouse Monoclonal antibody. Positive WB detected in human brain tissue. Observed molecular weight by Western-blot: 37 kDa

## Preparation

This antibody was obtained by immunization of IL-17F recombinant protein (Accession Number: NM\_052872). Purification method: Protein G purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:107385(IL17F antibody) at dilution of 1:1000