

IL1F10 antibody



Catalog Number: 107387

Product name

IL1F10 antibody

Immunogen

[Human IL1F10 Recombinant protein \(His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

IL1F10 Mouse Monoclonal antibody. Positive WB detected in Recombinant protein.

Preparation

This antibody was obtained by immunization of IL1F10 recombinant protein (Accession Number: NM_173161). Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



Recombinant protein were subjected to SDS PAGE followed by western blot with Catalog No:107387(IL1F10 Antibody) at dilution of 1:1000