

IL36RN antibody



Catalog Number: 107399

Product name

IL36RN antibody

Immunogen

Specificity

Human; other species not tested.

Antibody description

IL36RN Mouse Monoclonal antibody. Positive WB detected in Recombinant protein.

Preparation

This antibody was obtained by immunization of IL36RN recombinant protein (Accession Number: NM_012275). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2b

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



Recombinant protein were subjected to SDS PAGE followed by western blot with Catalog No:107399(IL1F5 Antibody) at dilution of 1:2000