

IL36RN antibody



Catalog Number: 107399

General Information

Product name:

IL36RN antibody

Immunogen:

Applications:

ELISA, WB

Specificity:

Human; other species not tested.

Antibody description:

IL36RN Mouse Monoclonal antibody. Positive WB detected in Recombinant protein.

Properties

Preparation:

This antibody was obtained by immunization of IL36RN recombinant protein (Accession Number: NM_012275). Purification method: Protein A purified.

Formulation:

PBS with 0.02% sodium azide and 50% glycerol

pH 7.3.

Storage:

Store at -20°C. DO NOT ALIQUOT

Clonality:

Monoclonal


Ig Type: Mouse IgG2b

Applications

Recommended Dilution:

WB: 1:1000-1:10000

Validations



Recombinant protein were subjected to SDS PAGE followed by western blot with Catalog No:107399(IL1F5 Antibody) at dilution of 1:2000