

Mcp1 antibody

Catalog Number: 107420

Product name

Mcp1 antibody

Specificity

Mouse; other species not tested.

Antibody description

Mcp1 Mouse Monoclonal antibody. Positive WB detected in LPS treated RAW 264.7 cells. Observed molecular weight by Western-blot: 16 kDa

Preparation

This antibody was obtained by immunization of Mcp1 recombinant protein (Accession Number: NM-011333). Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



LPS treated RAW 264.7 cells were subjected to SDS PAGE followed by western blot with Catalog No:107420(Mcp1 Antibody) at dilution of 1:1000