

# OTUD7B antibody



Catalog Number: 107439

## Product name

OTUD7B antibody

## Specificity

Human; other species not tested.

## Antibody description

OTUD7B Mouse Monoclonal antibody. Positive WB detected in Jurkat cells. Observed molecular weight by Western-blot: 92 kDa

## Preparation

This antibody was obtained by immunization of OTUD7B recombinant protein (Accession Number: BC020622). Purification method: Protein A purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2b

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:107439(OTUD7B Antibody) at dilution of 1:1000