

# p504S,AMACR antibody



Catalog Number: 107445

## Product name

p504S,AMACR antibody

## Specificity

human, mouse; other species not tested.

## Antibody description

p504S,AMACR Mouse Monoclonal antibody.  
Positive WB detected in HEK-293 cells. Observed molecular weight by Western-blot: 42kd

## Preparation

This antibody was obtained by immunization of p504S,AMACR recombinant protein (Accession Number: NM\_203382). Purification method: Protein G purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:107445(p504S,AMACR antibody) at dilution of 1:1000