

Pan Acetylation antibody

Catalog Number: 107451

Product name

Pan Acetylation antibody

Specificity

Human, Mouse, Rat, Dog; other species not tested.

Antibody description

Pan Acetylation Mouse Monoclonal antibody. Positive WB detected in HeLa cells. Positive IF detected in MDCK cells. Observed molecular weight by Western-blot: 17 kDa

Preparation

This antibody was obtained by immunization of Peptide. Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:20-1:200

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:107451(Pan Acetylation Antibody) at dilution of 1:1500



Immunofluorescent analysis of (10% Formaldehyde) fixed MDCK cells using (Pan Acetylation Antibody) at dilution of 1:150 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Mouse IgG(H+L)