

# PCMT1 antibody



Catalog Number: 107457

## Product name

PCMT1 antibody

## Immunogen

[Human PCMT1 Recombinant protein \(His tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

PCMT1 Mouse Monoclonal antibody. Positive IF detected in HeLa cells. Positive WB detected in Raji cells, HeLa cells. Observed molecular weight by Western-blot: 25-27 kDa

## Preparation

This antibody was obtained by immunization of PCMT1 recombinant protein (Accession Number: NM\_001360452). Purification method: Protein G purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:20-1:200

## Validations



Raji cells were subjected to SDS PAGE followed by western blot with Catalog No:107457(PCMT1 Antibody) at dilution of 1:1000



Immunofluorescent analysis of HeLa cells using Catalog No:107457(PCMT1 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Mouse IgG(H+L)