

# PDCD4 antibody

Catalog Number: 107460

## Product name

PDCD4 antibody

## Specificity

Human; other species not tested.

## Antibody description

PDCD4 Mouse Monoclonal antibody. Positive IHC detected in human breast cancer tissue, human gliomas tissue. Positive WB detected in HeLa cells, HEK-293 cells, Jurkat cells. Observed molecular weight by Western-blot: 56 kDa

## Preparation

This antibody was obtained by immunization of PDCD4 recombinant protein (Accession Number: NM\_014456). Purification method: Protein G purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications

ELISA, IHC, WB


## Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IHC: 1:20-1:200

## Validations



Immunohistochemistry of paraffin-embedded human breast cancer slide using Catalog No:107460(PDCD4 Antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:107460(PDCD4 Antibody) at dilution of 1:4000