

# PGM1 antibody

Catalog Number: 107467

## Product name

PGM1 antibody

## Specificity

Human; other species not tested.

## Antibody description

PGM1 Mouse Monoclonal antibody. Positive IF detected in HepG2 cells. Positive WB detected in Hepg2 cells, HEK-293 cells. Observed molecular weight by Western-blot: 61kDa

## Preparation

This antibody was obtained by immunization of PGM1 recombinant protein (Accession Number: NM\_002633). Purification method: Protein A purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2b

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:20-1:200

## Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:107467(PGM1 antibody) at dilution of 1:1000



Immunofluorescent analysis of HepG2 cells using Catalog No:107467(PGM1 Antibody) at dilution of 1:50 and Rhodamine-Goat anti-Mouse IgG