

# PR antibody

Catalog Number: 107481

## Product name

PR antibody

## Immunogen

[Human PR Recombinant protein \(His tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

PR Mouse Monoclonal antibody. Positive IHC detected in human breast cancer tissue.

## Preparation

This antibody was obtained by immunization of PR recombinant protein (Accession Number: NM\_000926). Purification method: Protein G purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications


ELISA, IHC

## Dilutions


Recommended Dilution:

IHC: 1:50-1:500

## Validations



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using (PR Antibody) at dilution of 1:3200 (under 10x lens). heat mediated antigen retrieved with Tris-EDTA buffer(pH9).



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using (PR Antibody) at dilution of 1:3200 (under 40x lens). heat mediated antigen retrieved with Tris-EDTA buffer(pH9).