

Prealbumin/transthyretin antibody



Catalog Number: 107483

Product name

Prealbumin/transthyretin antibody

Specificity

Human; other species not tested.

Antibody description

Prealbumin/transthyretin Mouse Monoclonal antibody. Positive IHC detected in human liver tissue. Positive WB detected in Human brain. Observed molecular weight by Western-blot: 16 kDa

Preparation

This antibody was obtained by immunization of Prealbumin/transthyretin recombinant protein (Accession Number: NM_000371). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2a

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:107483(TTR antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human liver tissue slide using Catalog No:107483(TTR Antibody) at dilution of 1:50