

# RACGAP1 antibody



Catalog Number: 107496

## Product name

RACGAP1 antibody

## Specificity

Human; other species not tested.

## Antibody description

RACGAP1 Mouse Monoclonal antibody. Positive WB detected in HeLa cells, HEK-293 cells, human brain tissue. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 71 kDa

## Preparation

This antibody was obtained by immunization of RACGAP1 recombinant protein (Accession Number: XM\_006719359). Purification method: Protein G purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:20-1:200

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:107496(RACGAP1 antibody) at dilution of



1:1000

Immunofluorescent analysis of HeLa cells, using RACGAP1 antibody Catalog No: at 1:50 dilution and Rhodamine-labeled goat anti-mouse IgG (red).