

# RBPJ antibody



Catalog Number: 107503

---

## Product name

RBPJ antibody

## Immunogen

[Human RBPJ Recombinant protein \(His tag\)](#)

## Specificity

human; other species not tested.

## Antibody description

RBPJ Mouse Monoclonal antibody. Positive WB detected in HepG2 cells. Observed molecular weight by Western-blot: 60 kDa

## Preparation

This antibody was obtained by immunization of RBPJ recombinant protein (Accession Number: XM\_017008174). Purification method: Protein G purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG1


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:500-1:5000

## Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:107503(RBPJ antibody) at dilution of 1:1000