

S100A10 antibody



Catalog Number: 107517

Product name

S100A10 antibody

Specificity

Human; other species not tested.

Antibody description

S100A10 Mouse Monoclonal antibody. Positive WB detected in human lung tissue. Observed molecular weight by Western-blot: 11 kDa

Preparation

This antibody was obtained by immunization of S100A10 recombinant protein (Accession Number: NM_002966). Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



human lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:107517(S100A10 Antibody) at dilution of 1:2000