

SDCCAG8 antibody



Catalog Number: 107521

Product name

SDCCAG8 antibody

Specificity

Human; other species not tested.

Antibody description

SDCCAG8 Mouse Monoclonal antibody. Positive WB detected in HEK-293 cells, COLO 320 cells, THP-1 cells. Observed molecular weight by Western-blot: 83 kDa

Preparation

This antibody was obtained by immunization of SDCCAG8 recombinant protein (Accession Number: XM_005273023). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2b

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:107521(SDCCAG8 Antibody) at dilution of 1:1000