

ZAP70 antibody

Catalog Number: 107547

Product name

ZAP70 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

ZAP70 Mouse Monoclonal antibody. Positive WB detected in Jurkat cells. Observed molecular weight by Western-blot: 70 kDa

Preparation

This antibody was obtained by immunization of ZAP70 recombinant protein (Accession Number: BC039039). Purification method: .

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:107547(ZAP70 antibody) at dilution of 1:1000