

# SFRP2 antibody



Catalog Number: 107585

## Product name

SFRP2 antibody

## Specificity

Human, Rat; other species not tested.

## Antibody description

SFRP2 Mouse Monoclonal antibody. Positive WB detected in rat brain tissue, COLO 320 cells, Y79 cells. Observed molecular weight by Western-blot: 38 kDa

## Preparation

This antibody was obtained by immunization of SFRP2 recombinant protein (Accession Number: NM\_003013). Purification method: Protein G purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

## Validations



rat brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:107585(SFRP2 Antibody) at dilution of 1:2000