

SRD5A1 antibody



Catalog Number: 107600

Product name

SRD5A1 antibody

Specificity

Human; other species not tested.

Antibody description

SRD5A1 Mouse Monoclonal antibody. Positive WB detected in human testis tissue. Observed molecular weight by Western-blot: 32 kDa

Preparation

This antibody was obtained by immunization of SRD5A1 recombinant protein (Accession Number: NM_001047). Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:107600(SRD5A1 Antibody) at dilution of 1:1000