

# SRPX2 antibody



Catalog Number: 107602

## Product name

SRPX2 antibody

## Specificity

Human, Pig; other species not tested.

## Antibody description

SRPX2 Mouse Monoclonal antibody. Positive WB detected in pig lung tissue, A549 cells. Observed molecular weight by Western-blot: 53 kDa

## Preparation

This antibody was obtained by immunization of SRPX2 recombinant protein (Accession Number: NM\_014467). Purification method: Protein A purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2b

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



pig lung tissue were subjected to SDS PAGE followed by western blot with (SRPX2 Antibody) at dilution of 1:1000