

TDP-43 (human specific) antibody



Catalog Number: 107618

Product name

TDP-43 (human specific) antibody

Specificity

Human; other species not tested.

Antibody description

TDP-43 (human specific) Mouse Monoclonal antibody. Positive IHC detected in human brain(FTLD-U) tissue, human brain tissue, human gliomas tissue, human pancreas cancer tissue, human pancreas tissue. Positive FC detected in MCF-7 cells. Positive WB detected in K-562 cells, human brain tissue, human heart tissue, human placenta tissue. Positive IP detected in K-562 cells. Observed molecular weight by Western-blot: 43 kDa

Preparation

This antibody was obtained by immunization of TDP-43 (human specific) recombinant protein (Accession Number: BC001487). Purification method: Caprylic acid/ammonium sulfate precipitation.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

ELISA, IP, FC, WB, IHC

Dilutions


Recommended Dilution:

WB: 1:5000-1:50000


IP: 1:500-1:5000

IHC: N/A

Validations





40X of FTLD-U case stained by Catalog No:115925 and Catalog No:107618, showing dystrophic neurites. (Figs were provided by Linda K. Kwong)



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:107618(TDP-43 Antibody) at dilution of 1:10000

TDP-43 (human specific) antibody

Catalog Number: 107618



IP Result of anti-TDP-43 (IP:Catalog No:107618, 5ug; Detection:Catalog No:107618 1:1000) with K-562 cells lysate 1720ug.

1×10^6 MCF-7 cells were stained with 0.2ug TDP-43 antibody (Catalog No:107618, red) and control antibody (blue). Fixed with 90% MeOH blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Mouse IgG(H+L) with dilution 1:1500.