

Thyroglobulin antibody

Catalog Number: 107623

Product name

Thyroglobulin antibody

Specificity

Human; other species not tested.

Antibody description

Thyroglobulin Mouse Monoclonal antibody. Positive WB detected in Transfected HEK-293 cells. Positive IHC detected in human thyroid tissue, human thyroid cancer tissue.

Preparation

This antibody was obtained by immunization of Thyroglobulin recombinant protein (Accession Number: NM_003235). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2a

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:50-1:500


Validations



Immunohistochemistry of paraffin-embedded human thyroid tissue slide using Catalog No:107623(TG Antibody) at dilution of 1:200 (under 10x lens).



Immunohistochemistry of paraffin-embedded human thyroid tissue slide using Catalog No:107623(TG Antibody) at dilution of 1:200 (under 40x lens).



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:107623(TG Antibody) at dilution of 1:1000