

# TPI1 antibody



Catalog Number: 107633

## Product name

TPI1 antibody

## Specificity

Human; other species not tested.

## Antibody description

TPI1 Mouse Monoclonal antibody. Positive WB detected in HEK-293 cells, HeLa cells. Observed molecular weight by Western-blot: 27kd

## Preparation

This antibody was obtained by immunization of TPI1 recombinant protein (Accession Number: NM\_000365). Purification method: Caprylic Acid/Ammonium Sulfate Precipitation.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:107633(TPI1 antibody) at dilution of 1:500