

TTL antibody

Catalog Number: 107650

Product name

TTL antibody

Specificity

Human; other species not tested.

Antibody description

TTL Mouse Monoclonal antibody. Positive WB detected in SH-SY5Y cells. Positive IHC detected in human liver. Observed molecular weight by Western-blot: 43kD

Preparation

This antibody was obtained by immunization of TTL recombinant protein (Accession Number: NM_153712). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2b

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:107650(TTL antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human liver using Catalog No:107650(TTL antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human liver using Catalog No:107650(TTL antibody) at dilution of 1:50 (under 40x lens)