

6*his,His tag antibody

Catalog Number: 107663

Product name

6*his,His tag antibody

Specificity

Tag

Antibody description

6*his,His tag Rabbit Polyclonal antibody. Positive WB detected in Recombinant protein.

Preparation

This antibody was obtained by immunization of Recombinant Protein. Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:5000

Validations



Western blot of 6*His-tagged fusion protein with anti-6*HIS tag (Catalog No:117296) at various dilutions.