

# ACD antibody



Catalog Number: 107733

## Product name

Polyclonal

ACD antibody

## Ig Type

## Specificity

Rabbit IgG

Human; other species not tested.

## Applications

## Antibody description

ELISA, WB

ACD Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells. Observed molecular weight by Western-blot: 57-60 kDa

## Dilutions

Recommended Dilution:


## Preparation

WB: 1:500-1:5000

This antibody was obtained by immunization of ACD recombinant protein (Accession Number: BC016904). Purification method: Antigen affinity purified.

## Validations

## Formulation



PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:107733(ACD Antibody) at dilution of 1:1500

## Clonality