

# ACK1 antibody



Catalog Number: 107740

## Product name

ACK1 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

ACK1 Rabbit Polyclonal antibody. Positive WB detected in SH-SY5Y cells, A431 cells, HeLa cells, HepG2 cells, mouse brain tissue, HT-1080 cells. Positive IHC detected in human breast cancer. Observed molecular weight by Western-blot: 70 kDa

## Preparation

This antibody was obtained by immunization of ACK1 recombinant protein (Accession Number: BC028164). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC

## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


## Validations



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:107740(ACK1 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:107740(ACK1 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:107740(ACK1 antibody) at dilution of 1:50 (under 40x lens)