

# ACOT6 antibody



Catalog Number: 107746

## Product name

ACOT6 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

ACOT6 Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue, HeLa cells. Positive IHC detected in human kidney tissue. Observed molecular weight by Western-blot: 23kd

## Preparation

This antibody was obtained by immunization of ACOT6 recombinant protein (Accession Number: NM\_001365789). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

## Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:107746(ACOT6 antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human kidney slide using Catalog No:107746(ACOT6 Antibody) at dilution of 1:50