

# ACVR1B antibody



Catalog Number: 107753

## Product name

ACVR1B antibody

## Immunogen

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

ACVR1B Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, DMSO treated HepG2 cells. Observed molecular weight by Western-blot: 49 kDa

## Preparation

This antibody was obtained by immunization of ACVR1B recombinant protein (Accession Number: NM\_004302). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:1000

## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:107753(ACVR1B Antibody) at dilution of 1:300