

ACVR1C antibody



Catalog Number: 107754

Product name

ACVR1C antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

ACVR1C Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, human brain tissue, human kidney tissue. Observed molecular weight by Western-blot: 55kd

Preparation

This antibody was obtained by immunization of ACVR1C recombinant protein (Accession Number: NM_145259). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:107754(ACVR1C antibody) at dilution of 1:300