

ADCY6 antibody

Catalog Number: 107781

Product name

ADCY6 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ADCY6 Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in human brain tissue, mouse brain tissue. Observed molecular weight by Western-blot: 131kd

Preparation

This antibody was obtained by immunization of ADCY6 recombinant protein (Accession Number: BC064923). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:1000

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:107781(ADCY6 antibody) at dilution of 1:200



IP Result of anti-ADCY6 (IP:Catalog No:107781, 4ug; Detection:Catalog No:107781 1:300) with mouse brain tissue lysate 4000ug.