

ADH1C antibody



Catalog Number: 107787

Product name

ADH1C antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ADH1C Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue, rat liver tissue. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 37-45 kDa

Preparation

This antibody was obtained by immunization of ADH1C recombinant protein (Accession Number: NM_000669). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IF: 1:10-1:100

Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:107787(ADH1C antibody) at dilution of 1:4000



Immunofluorescent analysis of HepG2 cells, using ADH1C antibody Catalog No:107787 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red). Blue pseudocolor = DAPI (fluorescent DNA dye).