

ADIPOR2 antibody



Catalog Number: 107794

Product name

ADIPOR2 antibody

Immunogen

[Human ADIPOR2 Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ADIPOR2 Rabbit Polyclonal antibody. Positive WB detected in L02 cells, HeLa cells, mouse uterus tissue. Positive IP detected in HeLa cells. Observed molecular weight by Western-blot: 44-50 kDa

Preparation

This antibody was obtained by immunization of ADIPOR2 recombinant protein (Accession Number: XM_047429546). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



L02 cells were subjected to SDS PAGE followed by western blot with Catalog No:107794(ADIPOR2 antibody) at dilution of 1:500



IP Result of anti-ADIPOR2 (IP:Catalog No:107794, 4ug; Detection:Catalog No:107794 1:500) with HeLa cells lysate 2800ug.