

ALKBH5 antibody



Catalog Number: 107858

Product name

ALKBH5 antibody

Immunogen

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ALKBH5 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HEK-293 cells, HepG2 cells. Observed molecular weight by Western-blot: 40-50 kDa

Preparation

This antibody was obtained by immunization of ALKBH5 recombinant protein (Accession Number: BC062339). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:1000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:107858(ALKBH5 antibody) at dilution of 1:400