

Alpha galactosidase A antibody



Catalog Number: 107871

Product name

Alpha galactosidase A antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Alpha galactosidase A Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells. Positive IF detected in HepG2. Observed molecular weight by Western-blot: 49kd

Preparation

This antibody was obtained by immunization of Alpha galactosidase A recombinant protein (Accession Number: NM_000169). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:107871 (GLA antibody) at dilution of 1:1500



Immunofluorescent analysis of HepG2 cells, using GLA antibody Catalog No:107871 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).