

# ALX3 antibody



Catalog Number: 107887

## Product name

ALX3 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

ALX3 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, HEK-293 cells. Observed molecular weight by Western-blot: 70 kDa

## Preparation

This antibody was obtained by immunization of ALX3 recombinant protein (Accession Number: NM\_006492). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:107887(ALX3 Antibody) at dilution of 1:1000