

ADRA1B-Specific antibody



Catalog Number: 107890

Product name

ADRA1B-Specific antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ADRA1B-Specific Rabbit Polyclonal antibody.
Positive IHC detected in mouse brain tissue,
human heart tissue.

Preparation

This antibody was obtained by immunization of
Peptide (Accession Number: NM_000679).
Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH
7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:50-1:500

Validations



Immunohistochemistry of paraffin-embedded
mouse brain tissue slide using Catalog
No:107890(ADRA1B-Specific Antibody) at dilution
of 1:100 (under 10x lens). heat mediated antigen
retrieved with Tris-EDTA buffer(pH9).



Immunohistochemistry of paraffin-embedded
mouse brain tissue slide using Catalog
No:107890(ADRA1B-Specific Antibody) at dilution
of 1:100 (under 40x lens). heat mediated antigen
retrieved with Tris-EDTA buffer(pH9).