

# AKAP11-Specific antibody



Catalog Number: 107929

## Product name

AKAP11-Specific antibody

## Specificity

Human; other species not tested.

## Antibody description

AKAP11-Specific Rabbit Polyclonal antibody. Positive IHC detected in human kidney tissue, human brain tissue. Positive WB detected in human placenta tissue. Observed molecular weight by Western-blot: 220 kDa

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_016248). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC

## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200


## Validations



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:107929(AKAP11-Specific antibody) at dilution of 1:100



Immunohistochemistry of paraffin-embedded human kidney tissue slide using Catalog No:107929(AKAP11-Specific Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human kidney tissue slide using Catalog No:107929(AKAP11-Specific Antibody) at dilution of 1:50 (under 40x lens)