

AKR7A3 antibody



Catalog Number: 107946

Product name

AKR7A3 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

AKR7A3 Rabbit Polyclonal antibody. Positive WB detected in human liver tissue, HepG2 cells, human brain tissue, NIH/3T3 cells. Observed molecular weight by Western-blot: 37kd and 55-60kd

Preparation

This antibody was obtained by immunization of AKR7A3 recombinant protein (Accession Number: NM_012067). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:107946(AKR7A3 antibody) at dilution of 1:400