

# ALDH5A1 antibody



Catalog Number: 107972

## Product name

ALDH5A1 antibody

IHC: 1:20-1:200


## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

ALDH5A1 Rabbit Polyclonal antibody. Positive IHC detected in human liver tissue. Positive WB detected in mouse pancreas tissue, human heart tissue, human skeletal muscle tissue, mouse brain tissue, mouse kidney tissue. Observed molecular weight by Western-blot: 51kd

## Validations



mouse pancreas tissue were subjected to SDS PAGE followed by western blot with Catalog No:107972(ALDH5A1 antibody) at dilution of 1:500

## Preparation

This antibody was obtained by immunization of ALDH5A1 recombinant protein (Accession Number: NM\_001080). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:

WB: 1:500-1:5000



Immunohistochemistry of paraffin-embedded human liver tissue slide using Catalog No:107972(ALDH5A1 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human liver tissue slide using Catalog

# ALDH5A1 antibody



Catalog Number: 107972

---

No:107972(ALDH5A1 Antibody) at dilution of 1:50 (under 40x lens)